CREDO[®] heal Stent

Patienteninformation

Patient information

acandis[®]

ENGINEERING STROKE SOLUTIONS

Acandis GmbH

Theodor-Fahrner-Str. 6

75177 Pforzheim/Germany

www.acandis.com

CE 0297

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1 Product description

The CREDO® heal Stent is a self-expanding nitinol stent with a fibrin/heparin coating. The stent has three platinum-iridium radiopaque markers at each end.

2 Indication

Used in combination with the NeuroSpeed® PTA Balloon Catheter, the CREDO® heal Stent is indicated to increase the vessel diameter of intracranial arteries (for example the ICA, MCA, VA/V4, BA) that are suitable for and accessible to the system (stent with PTA balloon catheter) in patients:

- with arteriosclerotic (probable according to clinical, anamnestic or laboratory chemical testing) vessel constrictions (high-grade intracranial stenoses > 70 %) with significant circulatory disturbance (haemodynamic infarct pattern and/or indication of restricted collaterals),
 - with recurrent symptoms in spite of antiaggregation therapy (therapy-refractory) after at least two events (stroke),
 - over 18 years old,
 - whose last stroke occurred at least 7 days before and in whom the impairment after stroke has a rating of 3 or less on the modified Rankin scale (mRS) at the time of treatment.
- or
- with evidence of intracranial large vessel occlusion (LVO),
 - over 18 years old,
 - with acute persistent vascular occlusion due to high-grade intracranial stenosis for which alternative therapy concepts are not considered or fail,
 - with suspected underlying stenosis of the occluded artery, which is suitable for PTA and stenting according to the treating physician's assessment,
 - with small to moderate infarct core prior to initiation of thrombectomy (CT ASPECTS 6-10, DWI lesion < 70 ml).

3 Information on magnetic resonance imaging (MRI)

Non-clinical tests have shown that the CREDO® heal Stent is suitable for MR. After implantation of an CREDO® heal Stent, patients can be safely scanned with a static magnetic field of 3-Tesla. The rise in temperature of the CREDO® heal Stent determined in non-clinical tests during MR imaging for 15 minutes with a 3-Tesla MR system was less than 0.5 °C.

The MR imaging quality may be affected if the CREDO® heal Stent is located in the area of interest. Optimisation of the imaging parameters is recommended. Please bring this information together with your implant card to all future MRI appointments.

4 Further information

For more information about this implant as well as updates on the information in this leaflet, please visit our website at www.acandis.com/patientimplantinfo.

The CREDO® heal Stent remains implanted for life. In case of questions, medical concerns, health symptoms, as well as operating instructions and necessary follow-up examinations, please consult your attending physician or healthcare professional. In case of an emergency, please notify emergency medical services.