ACCERO[®] heal Stent

Patienteninformation

Patient information

acandis[®]

ENGINEERING STROKE SOLUTIONS

Acandis GmbH

Theodor-Fahrner-Str. 6

75177 Pforzheim/Germany

www.acandis.com

CE 0297

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1 Product description

The ACCERO® heal Stent is a self-expanding nitinol stent with a fibrin/heparin coating. The stent has platinum-iridium radiopaque markers

2 Indication

The ACCERO® heal Stent is intended for use with embolisation materials in the treatment of intracranial aneurysms.

3 Information on magnetic resonance imaging (MRI)

Non-clinical tests have shown that the ACCERO® heal Stent is suitable for MR. After implantation of an ACCERO® heal Stent patients can be safely scanned with a static magnetic field of 3 Tesla.

The MR imaging quality may be affected if the ACCERO® heal Stent is located in the area of interest. Optimisation of the imaging parameters is recommended.

The patient with the implanted ACCERO® heal Stent can undergo an MRI scan under the following conditions:

- o Static magnetic field strength 1.5 and 3.0 Tesla
- o Maximum spatial gradient field 7.0 T/m at 3T, and 11.0 T/m at 1.5T
- o Maximum MR System reported, whole-body-averaged Specific Absorption Rate (SAR) for 15 minutes of scanning 1.4 W/kg at 1.5T produced an additional maximum temperature increase of 0.5 °C during 15 minutes of MRI; at 3-Tesla the rise was less than 0.2°C
- o In non-clinical testing, the image artefact caused by the device extends approximately 2 mm from the implanted device when imaged with a gradient echo sequence and a 3 Tesla MRI system
- o In non-clinical testing, the magnetically induced displacement force and magnetically induced torque were tested and no clinically significant displacement or torque was measured

4 Further information

For more information about this implant as well as updates on the information in this leaflet, please visit our website at www.acandis.com/patientimplantinfo.

The ACCERO® heal Stent remains implanted for life. In case of questions, medical concerns, health symptoms, as well as operating instructions and necessary follow-up examinations, please consult your attending physician or healthcare professional. In case of an emergency, please notify emergency medical services.