ACCERO[®] Rex Stent

Patienteninformation
Patient information

acandis[®]

ENGINEERING STROKE SOLUTIONS

Acandis GmbH

Theodor-Fahrner-Str. 6

75177 Pforzheim/Germany

www.acandis.com

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1 Product information

The ACCERO® Rex Stent is a self-expanding stent made of nitinol. In the centre, the stent has a radiopaque marker made of platinum-iridium.

2 Indication

The ACCERO® Rex Stent is intended for use with embolisation materials in the treatment of intracranial aneurysms.

3 Information on Magnetic Resonance Imaging (MRI)

Non-clinical tests have shown that the ACCERO® Rex Stent is partly suitable for MR. A patient can be safely scanned after implantation of the ACCERO® Rex Stent under the following conditions:

- Stent length up to 60 mm
- Static magnetic field of max. 3.0 tesla
- Static magnetic field gradient of max. 25 T/m
- Maximum specific absorption rate indicated by the MRI system, averaged over the whole body (SAR) of 4 W/kg for a 15-minute scan
- Controlled operating mode of the MRI system (Level 1)

In non-clinical tests, a temperature rise of the ACCERO® Rex Stent of less than 4 °C was detected in MR imaging at the above-mentioned conditions for 15 minutes. MR image quality may be impaired if the affected area is in the same range as the ACCERO® Rex Stent. Optimisation of the imaging parameters is recommended.

Please bring this information together with your implant card to all future MRI appointments.

4 Further information

ACCERO® Rex Stent is implanted for life. ACCERO® Rex Stent is a passive implant that does not require maintenance or monitoring. Nevertheless, it is necessary for you to follow the instructions of your treating physician. These include taking the prescribed medications and attending arranged follow-up examinations. In combination with a healthy lifestyle, this can reduce the risk of recurring symptoms.

If you have any questions, medical concerns or health-related symptoms and with regard to recommended action and required follow-up examinations, please consult your treating physician or your specialist care team.

In an emergency, please notify the emergency medical service.

Your treating physician will inform you when you should seek medical assistance. Furthermore, he will let you know of any potential side effects that might occur in connection with the implant, such as: unconsciousness, headaches, drowsiness, dizziness, speech disorders, uneven pupil dilation, nausea, fluctuations in blood pressure, mental confusion.

Please discuss with your treating physician which of these side effects might occur in your case.

This implant was manufactured in consideration of all requirements relating to medical devices. Biocompatibility tests have shown that any manufacturing residuals do not compromise the safety of the product.

5 Reporting of adverse events

Please consult your treating physician to report side effects and adverse events that might be related to the implant.

For more information on this implant and updates to the information on this printout, please visit our website at www.acandis.com/productinfo

You can report potential serious events in connection with this implant to Acandis®. Additionally, patients in Australia can contact the Australian Therapeutic Goods Administration (TGA) at: www.tga.gov.au