CARESTO[®] heal Stent

Patienteninformation
Patient information

acandis[®]

ENGINEERING STROKE SOLUTIONS

Acandis GmbH

Theodor-Fahrner-Str. 6

75177 Pforzheim/Germany

www.acandis.com

CE 0297

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1 Product Information

The CARESTO® heal Stent is a self-expanding stent made of nitinol (nickel-titanium alloy) with a platinum core and fibrin/heparin coating. The stent is fully visible under fluoroscopy.

2 Intended Use

2.1 Intended Purpose

The purpose of the CARESTO® heal Stent is to enlarge the vascular diameter of the common carotid artery (CCA), the internal carotid artery (ICA) and carotid bifurcation.

2.2 Indication

The CARESTO® heal Stent is indicated for treatment of stenoses in the common carotid artery (CCA), the internal carotid artery (ICA) and carotid bifurcation.

3 Information on Magnetic Resonance Imaging (MRI)

Non-clinical testing has demonstrated the CARESTO® heal Stent is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 and 3.0 T
- Maximum spatial field gradient of 11.0 T/m at 1.5 T and 7.0 T/m at 3.0 T
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg at 1.5 T and 2.0 W/kg at 3.0 T

Under the scan conditions defined above, the CARESTO® heal Stent is expected to produce a maximum temperature rise of less than 0.7 °C at 3.0 T and less than 1.6 °C at 1.5 T after 15 minutes of continuous scanning.

In non-clinical testing, the image artefact caused by the device extends up to 5 mm from the CARESTO® heal Stent when imaged with a spin echo pulse sequence and a 3.0 T MR system or up to 3 mm when imaged with a gradient echo pulse sequence and a 3.0 T MR system. Optimisation of the imaging parameters is recommended.

The MRI compatibility of the implant has only been established with a single stent. Acandis® offers no guarantee as to the MRI compatibility of overlapping stents.

Please bring this information together with your implant card to all future MRI appointments.

4 Further information

The CARESTO® heal Stent is implanted for life. The CARESTO® heal Stent is a passive implant that does not require maintenance or monitoring. Nevertheless, it is necessary for you to follow the instructions of your treating physician. These include taking the prescribed medications and attending arranged follow-up examinations. In combination with a healthy lifestyle, this can reduce the risk of recurring symptoms.

If you have any questions, medical concerns or health-related symptoms and with regard to recommended action and required follow-up examinations, please consult your treating physician or your specialist care team.

In an emergency, please notify the emergency medical service.

Your treating physician will inform you when you should seek medical assistance. Furthermore, he will let you know of any potential side effects that might occur in connection with the implant, such as: unconsciousness, headaches, drowsiness, dizziness, speech disorders, uneven pupil dilation, nausea, fluctuations in blood pressure, mental confusion.

Please discuss with your treating physician which of these side effects might occur in your case.

This implant was manufactured in consideration of all requirements relating to medical devices. Biocompatibility tests have shown that any manufacturing residuals do not compromise the safety of the product.

5 Reporting of adverse events

Please consult your treating physician to report side effects and adverse events that might be related to the implant.

For more information on this implant and updates to the information on this printout, please visit our website at www.acandis.com/productinfo

You can report potential serious events in connection with this implant to Acandis®. Additionally, patients in Australia can contact the Australian Therapeutic Goods Administration (TGA) at: www.tga.gov.au