CREDO® Stent

DE	Patienteninformation
EN	Patient information
FR	Informations relatives au patient
NL	Patiënteninformatie
IT	Informazioni ai pazienti
ES	Información para el paciente
PT	Informação ao paciente
FI	Potilastiedot
SV	Patientinformation
DA	Patientinformation
NO	Pasientinformasjon
TR	Hasta bilgilendirme broşürü
ET	Patsienditeave
LV	Informācija par pacientu
LT	Informācija par pacientu
PL	Informacja dla pacjenta
CS	Informace pro pacienta
SK	Informácie pre pacienta
SL	Informacije o bolniku
HR	Informacija za pacijenta
HU	Betegtájékoztató
RO	Informații pentru pacient
BG	Информация за пациента
EL	Πληροφορίες για τον ασθενή

ENGINEERING STROKE SOLUTIONS

Acandis GmbH

Theodor-Fahrner-Str. 6
75177 Pforzheim/Germany
www.acandis.com

	Deutsch	English	Français	Nederlands	Italiano	Español	Português	Suomi	Svensk	Dansk	Norsk	Türkçe
	Patienten-identifikation	Patient identification	Identification du patient	Patiënt-ID	Identificazione del paziente	Identificación del paciente	Identificação do paciente	Potilastunniste	Patientidentifiering	Patient-identifikation	Pasient-identifikasjon	Hasta kimliği
	Implantations-datum	Implantation date	Date d'implantation	Implantatiedatum	Data di impianto	Fecha de implante	Data do implante	Implantoinnin päivämäärä	Implantations-datum	Implantationsdato	Implantasjonsdato	İmplantasyon tarihi
	Name und Adresse der Gesundheits-einrichtung	Name and address of the healthcare facility	Nom et adresse de l'établissement de santé	Naam en adres van de gezondheidsz-orginstelling	Nome e indirizzo della struttura sanitaria	Nombre y dirección del centro sanitario	Nome e morada da instituição de saúde	Terveystietojen yksikön nimi ja osoite	Namn och adress för vårdinrättningen	Sundheds-institutionens navn og adresse	Helseinstitusjonens navn og adresse	Sağlık kurumunun unvanı ve adresi
	Hersteller	Manufacturer	Fabricant	Fabrikant	Fabbricante	Fabricante	Fabricante	Valmistaja	Tillverkare	Producent	Produsent	Üretici
	Internetseite mit Patienten-informationen	Website with patient information	Site Internet avec informations destinées aux patients	Internetpagina met informatie voor de patiënt	Sito web con informazioni per i pazienti	Sito web con información para el paciente	Página de internet com informações ao paciente	Potilastietosivusto	Internetsajt med patientinformation	Internetside med patientinformation	Nettside med pasientinformasjon	Hasta bilgilendirme broşürünün bulunduğu internet sitesi
	Medizinprodukt	Medical Device	Dispositif médical	Medisch hulpmiddel	Dispositivo medico	Producto médico	Dispositivo médico	Lääkinnällinen laite	Medicinteknisk produkt	Medicinsk produkt	Medisinsk utstyr	Tıbbi ürün
	Katalognummer	Catalogue number	Référence catalogue	Catalogusnummer	Numero di catalogo	Número de catálogo	Número do catálogo	Luettelonumero	Katalognummer	Katalognummer	Katalognummer	Katalog numarası
	Chargen-bezeichnung	Batch code	LOT	Batchbenaming	Numero di lotto	Denominación de lote	Designação do lote	Erämerkintä	Batchkod	Batch-betegnelse	Partikode	Parti tanımı
	Eindeutige Produkt-identifizierung	Unique Device Identifier	Identifiant unique du dispositif	Unieke product-identificatie	Identificazione univoca dei dispositivi	Identificación única del producto	Identificação clara do produto	Selkeä tuotetunniste	Unik produktiden-tifiering	Unik produkt-identifikation	Unik utstyrs-identifikasjon	Benzersiz ürün tanımlaması
	MR bedingt	MR Conditional	RM conditionnelle	MR voorwaardelijk	Compatibilità magnetica condizionata	Condicional para RM	Condicionado para RM	Tietyin edellytyksin MR-turvallinen	MR-villikorig	MR betinget	MR-betinget	Belirli koşullarda MR uyumlu

	Eesti	Latviešu	Lietuvių	Polski	Čeština	Slovenčina	Slovenščina	Hrvatski	Magyar	Română	Български	Ελληνικά
	Patsiendi tuvastamine	Pacienta identifikācija	Pacientų tapatybės nustatymas	Identyfikacja pacjenta	Identifikace pacienta	Identifikator pacienta	Identifikacija bolnika	Identifikacija pacijenta	Betegazonosító	Identificarea pacientului	Идентификация на пациента	Τυποποίηση ασθενούς
	Implanteerimise kuupäev	Implantācijas datums	Implantacijos data	Data implantacji	Datum implantace	Dátum implantácie	Datum vsaditve	Datum ugradnje	Implantáció dátuma	Data implantării	Дата на имплантиране	Ημερομηνία εμφύτευσης
	Tervishoiuasutuse nimi ja aadress	Veselības iestādes nosaukums un adrese	Sveikatos priežiūros įstaigos pavadinimas ir adresas	Nazwisko i adres placówki medycznej	Název a adresa zdravotnického zařízení	Názov a adresa zdravotníckeho zariadenia	Ime in naslov zdravstvene ustanove	Naziv i adresa zdravstvene ustanove	Egészségügyi intézmény neve és címe	Numele și adresa unității medicale	Наименование и адрес на здравното заведение	Όνομα και διεύθυνση της μονάδας υγείας
	Tootja	Ražotājs	Gamintojas	Producent	Výrobce	Proizvajalec	Výrobca	Proizvođač	Gyártó	Producător	Производител	Κατασκευαστής
	Veebisait patsiendi teabega	Interneta vietne ar informāciju par pacientu	Interneto svetainė su informacija apie pacientus	Strona internetowa z informacjami dla pacjenta	Internetová stránka s informacemi o pacientovi	Internetová stránka s informáciami pre pacienta	Spletna stran z informacijami za bolnike	Internetska stranica s informacijama za pacijenta	Weboldal betegtájékoztatókkal	Pagina de internet cu informații pentru pacient	Интернет страница с информация за пациента	Ιστότοπος με πληροφορίες για τον ασθενή
	Meditsiini toode	Medicīniskis izstrādājums	Medicininis gaminyis	Wyrób medyczny	Lékařský produkt	Medicinski pripomoček	Zdravotniška pomočka	Medicinski proizvod	Orvostechnikai eszköz	Produs medicinal	Медицинско изделие	Ιατροτεχνολογικό προϊόν
	Katalooginumber	Kataloga numurs	Katalogo numeris	Numer katalogowy	Číslo výrobku	Kataloška številka	Číslo položky	Kataloški broj	Cikkszám	Numărul articolului	Артикулен номер	Αριθμός προϊόντος
	Partii tähis	Partijas nosaukums	Partijos pavadinimas	Nazwa partii	Kód šarže	Poimenovanje šarže	Kód šarže	Broj serije	Tételkód	Cod de lot	Код на партидата	Κωδικός παρτίδας
	Toote ühene identifitseerimine	Skaidra izstrādājuma identifikācija	Unikalus gaminio identifikatorius	Jednoznaczna identyfikacja wyrobu	Jednoznačná identifikace produktu	Edinstvena identifikacija izdelka	Unikatny identifikator pomočky	Jedinstvena identifikacija proizvoda	Egyedi termékazonosító	Identificarea clară a produsului	Унікален ідентифікатор на izdelieto	Μοναδικός αναγνωριστικός αριθμός προϊόντος
	Osaliselt MR tundlik	MR nosacīti	Su išlygomis galimas MR	Warunkowo bezpieczny w środowisku MR	Podmíněno MR	MR pogojeno	Podmienneče bezpečné v prostredí MR	Uvjetno sigurno prilikom snimanja MR-om	MRI közelében korlátozottan használható	Condiționat RM	Условно съвместим с MP	Για ΜΣ

DE	Patienteninformation	5
EN	Patient information	7
FR	Informations relatives au patient	9
NL	Patiënteninformatie	11
IT	Informazioni ai pazienti	13
ES	Información para el paciente	15
PT	Informação ao paciente	17
FI	Potilastiedot	19
SV	Patientinformation	21
DA	Patientinformation	23
NO	Pasientinformasjon	25
TR	Hasta bilgilendirme broşürü	27
ET	Patsienditeave	29
LV	Informācija par pacientu	31
LT	Informacija apie pacientą	33
PL	Informacja dla pacjenta	35
CS	Informace pro pacienta	37
SK	Informácie pre pacienta	39
SL	Informacije o bolniku	41
HR	Informacija za pacijenta	42
HU	Betegtájékoztató	44
RO	Informații pentru pacient	46
BG	Информация за пациента	48
EL	Πληροφορίες για τον ασθενή	50

1 Produktinformation

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1 Product Information

The CREDO® Stent is a self-expanding stent made of nitinol. The stent has three radiopaque markers made of platinum-iridium at each end.

Patients may come into contact with the following implant materials:

Nickel-titanium alloy – strut width and height ≤ 0.1 mm,
max. 12 struts in cross-section

Platinum-iridium marker – 6 markers

2 Intended Use

2.1 Intended purpose

The CREDO® Stent is intended for the treatment of intracranial stenosis by keeping the blood vessel open after necessary pretreatment has been performed.

2.2 Indication

The CREDO® Stent is indicated, together with the NeuroSpeed® PTA Balloon Catheter, for the dilatation of the diameter of intracranial arteries (for example ICA, MCA, VA/V4 or BA) which are accessible to and suitable for the system (stent with PTA balloon catheter), in patients:

- with stenosis (high-grade intracranial stenosis of $> 70\%$) probably resulting from arteriosclerosis (as established clinically, from case history or from laboratory tests), where there is clear circulatory disturbance (haemodynamic infarction pattern and evidence of limited collateral vessels) despite intensive drug therapy,
- who have become symptomatic again despite antiplatelet therapy following at least two events (stroke) (treatment-refractory),
- who are over 18 years of age,
- in whom the most recent stroke occurred at least 7 days previously, with a diagnosed impairment after the stroke of 3 or less on the modified Rankin Scale (mRS) at the time of treatment.

or

- with evidence of an intracranial large vessel occlusion (LVO),
- who are over 18 years of age,
- with an acute, persistent vessel occlusion resulting from a high-grade intracranial stenosis, for which alternative treatment options are not possible or have failed,
- with a suspected underlying stenosis of the occluded artery, which is amenable to PTA and stenting in the opinion of the treating physician,
- with a small to medium infarct core prior to the commencement of thrombectomy (CT score ASPECTS 6-10, DWI lesion < 70 ml).

3 Information on Magnetic Resonance Imaging (MRI)

Non-clinical testing has demonstrated the CREDO® Stent is MR Conditional. A patient with this implant can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial magnetic field gradient of 4.2 T/m at 1.5 T and 11.0 T/m at 3.0 T
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.9 W/kg (at 1.5 T)/3.7 W/kg (at 3.0 T)

Under the scan conditions defined above, the CREDO® Stent is expected to produce a maximum temperature rise of less than 0.5°C (at 1.5 T)/ 0.5°C (at 3.0 T) after 15 minutes of continuous scanning.

In non-clinical testing, the image artefact caused by the implant extends up to 3 mm from the CREDO® Stent when imaged with a gradient echo or spin echo pulse sequence and a 3 T MR system. Optimisation of the imaging parameters is recommended.

The MR compatibility of the implant has only been established with a single stent. Acandis® offers no guarantee as to the MR compatibility of overlapping stents.

Please bring this information together with your implant card to all future MRI appointments.

4 Further information

The CREDO® Stent is implanted for life. The CREDO® Stent is a passive implant that does not require maintenance or monitoring. Nevertheless, it is necessary for you to follow the instructions of your treating physician. These include taking the prescribed medications and attending arranged follow-up examinations. In combination with a healthy lifestyle, this can reduce the risk of recurring symptoms.

If you have any questions, medical concerns or health-related symptoms and with regard to recommended action and required follow-up examinations, please consult your treating physician or your specialist care team.

In an emergency, please notify the emergency medical service.

Your treating physician will inform you when you should seek medical assistance. Furthermore, he will let you know of any potential side effects that might occur in connection with the implant, such as: unconsciousness, headaches, drowsiness, dizziness, speech disorders, uneven pupil dilation, nausea, fluctuations in blood pressure, mental confusion.

Please discuss with your treating physician which of these side effects might occur in your case.

This implant was manufactured in consideration of all requirements relating to medical devices. Biocompatibility tests have shown that any manufacturing residuals do not compromise the safety of the product.

5 Reporting of adverse events

Please consult your treating physician to report side effects and adverse events that might be related to the implant.

For more information on this implant and updates to the information on this printout, please visit our website at: www.acandis.com/productinfo

You can report potential serious events in connection with this implant to Acandis®. Additionally, patients in Australia can contact the Australian Therapeutic Goods Administration (TGA) at: www.tga.gov.au

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