



BRILLIANT HIGH-FLYER

ACCERO® Stent

SIZING SUPPORT CHART

Labelled ACCERO® Dimensions (mm)	Reference Number		Unconstrained ACCERO® Dimensions (mm)	ACCERO® Lengths in corresponding Intended Use Diameters (mm)		
		Ø	2.7	2.5	2.0	1.5
2.5 x 10	01-000800	Stent Length	6	10	12	14
2.5 x 15	01-000801		10	15	18	22
2.5 x 20	01-000802		13	20	25	29
		Ø	3.2	3.0	2.5	2.0
3.0 x 10	01-000803	Stent Length	7	10	13	15
3.0 x 15	01-000804		10	15	19	22
3.0 x 20	01-000805		14	20	25	29
		Ø	3.7	3.5	3.0	2.5
3.5 x 10	01-000806	Stent Length	7	10	12	14
3.5 x 15	01-000807		10	15	19	22
3.5 x 20	01-000808		14	20	26	30
3.5 x 25	01-000841		17	25	32	37
		Ø	4.2	4.0	3.5	3.0
4.0 x 15	01-000810	Stent Length	10	15	17	18
4.0 x 20	01-000811		14	20	24	27
4.0 x 25	01-000845		17	23	27	30
		Ø	4.7	4.5	4.0	3.5
4.5 x 15	01-000813	Stent Length	11	15	20	22
4.5 x 20	01-000814		15	20	26	29
4.5 x 25	01-000842		20	25	30	34

Note: all indicated lengths can vary within a tolerance range of +/- 1mm

For optimal case preparation, Acandis also offers software-based 3D Sizing Support.

For further information please contact: clinical-support@acandis.com