



SIZING SUPPORT CHART | ACCERO®

Labelled ACCERO® Dimensions (mm)		Unconstrained ACCERO® Dimensions (mm)	ACCERO® Lengths in corresponding Intended Use Diameters (mm)			
			Ø	2.5	2.0	1.5
2.5 × 10 2.5 × 15 2.5 × 20	Ø	2.7	2.5	2.0	1.5	
	Stent Length	6	10	12	14	
		10	15	18	22	
		13	20	25	29	
3.0 × 10 3.0 × 15 3.0 × 20	Ø	3.2	3.0	2.5	2.0	
	Stent Length	7	10	13	15	
		10	15	19	22	
		14	20	25	29	
3.5 × 10 3.5 × 15 3.5 × 20 3.5 × 25	Ø	3.7	3.5	3.0	2.5	
	Stent Length	7	10	12	14	
		10	15	19	22	
		14	20	26	30	
4.0 × 15 4.0 × 20 4.0 × 25	Ø	4.2	4.0	3.5	3.0	
	Stent Length	10	15	17	18	
		14	20	24	27	
		17	25	27	30	

Note: all indicated lengths can vary within a tolerance range of +/- 1 mm

For optimal case preparation, Acandis® also offers software-based 3D Sizing Support.
For further information please contact: clinical-support@acandis.com