


FLEXIBLE.
SMOOTH.
SECURE.



ACCLINO® flex plus Stent

NEXT GENERATION

Self-expanding nitinol stent with enhanced radiopacity and j-shaped transport wire tip for comprehensive treatment options.

ORDERING INFORMATION

Labelled ACCLINO® flex plus Dimensions (mm)	Reference Number	Stent Diameter (mm)	Stent Length (mm)	Recommended Vessel Diameter (mm)	Required Microcatheter for Delivery** (inch)
3.0 × 15	01-001122	3.0	15	1.5 – 2.5	0.0165 – 0.017
3.0 × 20	01-001123		20		
3.0 × 25	01-001124		25		
3.0 × 30	01-001125		30		
3.0 × 35	01-001126		35		
3.5 × 15	01-001132	3.5	15	1.5 – 3.0	0.0165 – 0.017
3.5 × 20	01-001133		20		
3.5 × 25	01-001134		25		
3.5 × 30	01-001135		30		
3.5 × 35	01-001136		35		
4.0 × 15*	01-001142	4.0	15	2.5 – 3.5	0.0165 – 0.017
4.0 × 20*	01-001143		20		
4.0 × 25*	01-001144		25		
4.0 × 30*	01-001145		30		
4.0 × 35*	01-001146		35		
4.5 × 15	01-001152	4.5	15	2.5 – 4.0	0.0165 – 0.017
4.5 × 20	01-001153		20		
4.5 × 25	01-001154		25		
4.5 × 30	01-001155		30		
4.5 × 35	01-001156		35		
5.5 × 20	01-001173	5.5	20	3.5 – 5.0	0.0165 – 0.017
5.5 × 25	01-001174		25		
5.5 × 30	01-001175		30		
5.5 × 35	01-001176		35		
6.5 × 20	01-001193	6.5	20	4.0 – 6.0	0.021
6.5 × 25	01-001194		25		
6.5 × 30	01-001195		30		
6.5 × 35	01-001196		35		

HRF: High Radial Force – compared to stents within the same recommended vessel diameter.

* Coming soon. For availability please contact your local representative from Acandis®.

** The ACCLINO® flex plus Stent is compatible with the Acandis® NeuroSlider® Microcatheter.

All changes or modifications, may be technical or other, or changes in the availability of products are expressly reserved.

Distributed by:

acandis®
ENGINEERING STROKE SOLUTIONS

ACANDIS GmbH
Theodor-Fahrner-Str. 6
75177 Pforzheim
Germany

Tel: +49 7231 155 00 0
Fax: +49 7231 155 00 129
E-Mail: info@acandis.com
www.acandis.com

CE 0297