



MINI BUT GREAT

DERIVO® mini
Embolisation Device

SIZING SUPPORT CHART

Labelled DERIVO® mini Dimensions (mm)	Reference Number	Unconstrained DERIVO® mini Dimensions (mm)	DERIVO® mini Lengths in corresponding Intended Use Diameters (mm)		
		Ø 2.7	2.5	2.0	1.5
2.5 x 15	01-000428	11	15	20	23
2.5 x 20	01-000429	14	21	29	34
2.5 x 20 w/o tip	01-000432	14	21	29	34
2.5 x 25	01-000430	19	27	38	45
2.5 x 25 w/o tip	01-000433	19	27	38	45
		Ø 3.2	3.0	2.5	2.0
3.0 x 15	01-000422	11	14	18	21
3.0 x 20	01-000423	15	20	27	31
3.0 x 20 w/o tip	01-000426	15	20	27	31
3.0 x 25	01-000424	19	26	36	42
3.0 x 25 w/o tip	01-000427	19	26	36	42
		Ø 3.7	3.5	3.0	2.5
3.5 x 15	01-000416	11	14	17	20
3.5 x 20	01-000417	14	19	24	28
3.5 x 20 w/o tip	01-000420	14	19	24	28
3.5 x 25	01-000418	19	24	32	37
3.5 x 25 w/o tip	01-000421	19	24	32	37

Note: all indicated lengths can vary within a tolerance range of +/- 1mm

For optimal case preparation, Acandis also offers software-based 3D Sizing Support.

For further information please contact: clinical-support@acandis.com



ENGINEERING STROKE SOLUTIONS