

SIZING SUPPORT CHART – DERIVO® EMBOLISATION DEVICE

Labelled DERIVO® Dimensions (mm)	Reference Number		Unconstrained DERIVO® Dimensions (mm)	DERIVO® Lengths in corresponding Intended Use Diameters (mm)		
		Ø	3.7	3.5	3.0	2.5
3.5 × 15	01-000408	Device Length	10	15	20	25
3.5 × 20	01-000409		13	20	27	32
3.5 × 25	01-000410		16	25	35	41
3.5 × 30	01-000411		19	30	41	48
3.5 × 40	01-000415		25	40	53	66
		Ø	4.2	4.0	3.5	3.0
4.0 × 15	01-000381	Device Length	11	15	20	25
4.0 × 20	01-000330		14	20	27	32
4.0 × 25	01-000335		17	25	35	41
4.0 × 30	01-000340		20	30	41	48
4.0 × 40	01-000360		26	40	53	66
		Ø	4.7	4.5	4.0	3.5
4.5 × 15	01-000382	Device Length	11	15	20	25
4.5 × 20	01-000331		14	20	27	32
4.5 × 25	01-000336		17	25	35	41
4.5 × 30	01-000341		20	30	41	48
4.5 × 40	01-000361		26	40	53	66
		Ø	5.2	5.0	4.5	4.0
5.0 × 15	01-000383	Device Length	11	15	20	23
5.0 × 20	01-000332		14	20	27	32
5.0 × 25	01-000337		17	25	35	41
5.0 × 30	01-000342		20	30	41	48
5.0 × 40	01-000362		26	40	53	62
5.0 × 50	01-000363		34	50	68	82
		Ø	5.7	5.5	5.0	4.5
5.5 × 15	01-000384	Device Length	11	15	20	23
5.5 × 20	01-000333		14	20	27	32
5.5 × 25	01-000338		17	25	35	41
5.5 × 30	01-000343		20	30	41	48
5.5 × 40	01-000364		26	40	53	62
5.5 × 50	01-000365		34	50	68	82
		Ø	6.2	6.0	5.5	5.0
6.0 × 15	01-000385	Device Length	11	15	20	23
6.0 × 20	01-000334		14	20	27	32
6.0 × 25	01-000339		17	25	35	41
6.0 × 30	01-000344		20	30	41	48
6.0 × 40	01-000366		26	40	53	62
6.0 × 50	01-000367		34	50	68	82

Note: all indicated lengths can vary within a tolerance range of +/- 1mm

For optimal case preparation, Acandis also offers software-based 3D Sizing Support.

For further information please contact the Clinical Support Team: clinical-support@acandis.com