



VISIBLE ADAPTABILITY

DERIVO®
Embolisation Device

SIZING SUPPORT CHART

Labelled DERIVO® Dimensions (mm)	Reference Number	Unconstrained DERIVO® Dimensions (mm)	DERIVO® Lengths in corresponding Intended Use Diameters (mm)			
		Ø 3.7	3.5	3.0	2.5	
3.5 × 15	01-000408	10	15	20	25	
3.5 × 20	01-000409	13	20	27	32	
3.5 × 25	01-000410	16	25	35	41	
3.5 × 30	01-000411	19	30	41	48	
3.5 × 40	01-000415	25	40	53	66	
		Ø 4.2	4.0	3.5	3.0	
4.0 × 15	01-000381	11	15	20	25	
4.0 × 20	01-000330	14	20	27	32	
4.0 × 25	01-000335	17	25	35	41	
4.0 × 30	01-000340	20	30	41	48	
4.0 × 40	01-000360	26	40	53	66	
		Ø 4.7	4.5	4.0	3.5	
4.5 × 15	01-000382	11	15	20	25	
4.5 × 20	01-000331	14	20	27	32	
4.5 × 25	01-000336	17	25	35	41	
4.5 × 30	01-000341	20	30	41	48	
4.5 × 40	01-000361	26	40	53	66	
		Ø 5.2	5.0	4.5	4.0	
5.0 × 15	01-000383	11	15	20	23	
5.0 × 20	01-000332	14	20	27	32	
5.0 × 25	01-000337	17	25	35	41	
5.0 × 30	01-000342	20	30	41	48	
5.0 × 40	01-000362	26	40	53	62	
5.0 × 50	01-000363	34	50	68	82	
		Ø 5.7	5.5	5.0	4.5	
5.5 × 15	01-000384	11	15	20	23	
5.5 × 20	01-000333	14	20	27	32	
5.5 × 25	01-000338	17	25	35	41	
5.5 × 30	01-000343	20	30	41	48	
5.5 × 40	01-000364	26	40	53	62	
5.5 × 50	01-000365	34	50	68	82	
		Ø 6.2	6.0	5.5	5.0	
6.0 × 15	01-000385	11	15	20	23	
6.0 × 20	01-000334	14	20	27	32	
6.0 × 25	01-000339	17	25	35	41	
6.0 × 30	01-000344	20	30	41	48	
6.0 × 40	01-000366	26	40	53	62	
6.0 × 50	01-000367	34	50	68	82	

Note: all indicated lengths can vary within a tolerance range of +/- 1mm

For optimal case preparation, Acandis also offers software-based 3D Sizing Support.

For further information please contact: clinical-support@acandis.com | www.acandis.com



ENGINEERING STROKE SOLUTIONS